

## Antibodies Product Datasheet

<b>Product Name:</b>	Bcl-x Antibody
<b>Catalog Number:</b>	B7027
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IF/ICC, WB, IHC-p, ELISA
<b>Dilutions:</b>	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human BCL-XL. AA range:28-77
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	BCL2L1
<b>Protein Name:</b>	Bcl-2-like protein 1
<b>Other Name:</b>	BCL2L1; BCL2L; BCLX; Bcl-2-like protein 1; Bcl2-L-1; Apoptosis regulator Bcl-X
<b>Human Swiss Prot No.:</b>	Q07817
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Bcl-x Polyclonal Antibody detects endogenous levels of Bcl-x protein.
<b>Storage and Stability:</b>	-20°C/1 year