

Antibodies Product Datasheet

Product Name:	SSTR1 Antibody
Catalog Number:	G744
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SSTR1. AA range:9-58
Concentration:	1 mg/ml
Gene Name:	SSTR1
Protein Name:	Somatostatin receptor type 1
Other Name:	SSTR1; Somatostatin receptor type 1; SS-1-R; SS1-R; SS1R; SRIF-2
Human Swiss Prot No.:	P30872
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SSTR1 Polyclonal Antibody detects endogenous levels of SSTR1 protein.
Storage and Stability:	-20°C/1 year