
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | CTCF Monoclonal Antibody |
| Catalog Number: | N0172 |
| Reactivity: | Human,Monkey |
| Applications: | WB,ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | Affinity purification |
| Immunogen: | Purified recombinant fragment of human CTCF expressed in E. Coli. |
| Concentration: | |
| Gene Name: | CTCF |
| Protein Name: | Transcriptional repressor CTCF |
| Other Name: | CTCF; Transcriptional repressor CTCF; 11-zinc finger protein; CCCTC-binding factor; CTCFL paralog |
| Human Swiss Prot No.: | P49711 |
| Formulation: | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol. |
| Specificity: | CTCF Monoclonal Antibody detects endogenous levels of CTCF protein. |
| Storage and Stability: | -20°C/1 year |