

Antibodies Product Datasheet

Product Name:	VEGF-D Antibody
Catalog Number:	R12-2408
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VEGF-D. AA range:153-202
Concentration:	1 mg/ml
Gene Name:	FIGF
Protein Name:	Vascular endothelial growth factor D
Other Name:	FIGF; VEGFD; Vascular endothelial growth factor D; VEGF-D; c-Fos-induced growth factor; FIGF
Human Swiss Prot No.:	O43915
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VEGF-D Polyclonal Antibody detects endogenous levels of VEGF-D protein.
Storage and Stability:	-20°C/1 year