

## Antibodies Product Datasheet

<b>Product Name:</b>	Rsk-1 Antibody
<b>Catalog Number:</b>	B0537
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human p90 RSK. AA range:539-588
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	RPS6KA1
<b>Protein Name:</b>	Ribosomal protein S6 kinase alpha-1
<b>Other Name:</b>	RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6 kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K; MAP kinase-activated protein kinase 1a; MAPK-activated protein kinase 1a; MAPKAP kinase 1a; MAPKAP
<b>Human Swiss Prot No.:</b>	Q15418
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Rsk-1 Polyclonal Antibody detects endogenous levels of Rsk-1 protein.
<b>Storage and Stability:</b>	-20°C/1 year