

## Antibodies Product Datasheet

<b>Product Name:</b>	DNAJA1 Antibody
<b>Catalog Number:</b>	C48352
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human DNAJA1
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	DNAJA1 DNAJ2 HDJ2 HSJ2 HSPF4
<b>Protein Name:</b>	DnaJ homolog subfamily A member 1 (DnaJ protein homolog 2) (HSDJ) (Heat shock 40 kDa protein 4) (Heat shock protein J2) (HSJ-2) (Human DnaJ protein 2) (hDj-2)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P31689
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of DNAJA1 at Human, Mouse,Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year