

Antibodies Product Datasheet

Product Name:	Claudin-3 Antibody
Catalog Number:	C0144
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Claudin 3. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	CLDN3
Protein Name:	Claudin-3
Other Name:	CLDN3; C7orf1; CPETR2; Claudin-3; Clostridium perfringens enterotoxin receptor 2; CPE-R 2; CPE-receptor 2; Rat ventral prostate.1 protein homolog; hRVP1
Human Swiss Prot No.:	O15551
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Claudin-3 Polyclonal Antibody detects endogenous levels of Claudin-3 protein.
Storage and Stability:	-20°C/1 year