

Antibodies Product Datasheet

Product Name:	NF-M Antibody
Catalog Number:	R12-2263
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NF-M. AA range:542-591
Concentration:	1 mg/ml
Gene Name:	NEFM
Protein Name:	Neurofilament medium polypeptide
Other Name:	NEFM; NEF3; NFM; Neurofilament medium polypeptide; NF-M; 160 kDa neurofilament protein; Neurofilament 3; Neurofilament triplet M protein
Human Swiss Prot No.:	P07197
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NF-M Polyclonal Antibody detects endogenous levels of NF-M protein.
Storage and Stability:	-20°C/1 year