

## Antibodies Product Datasheet

<b>Product Name:</b>	ApoA-V Monoclonal Antibody
<b>Catalog Number:</b>	N0033
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of human ApoA-V expressed in E. Coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	APOA5
<b>Protein Name:</b>	Apolipoprotein A-V
<b>Other Name:</b>	APOA5; RAP3; Apolipoprotein A-V; Apo-AV; ApoA-V; Apolipoprotein A5; Regeneration-associated protein 3
<b>Human Swiss Prot No.:</b>	Q6Q788
<b>Formulation:</b>	Purified antibody in PBS containing 0.03% sodium azide.
<b>Specificity:</b>	ApoA-V Monoclonal Antibody detects endogenous levels of ApoA-V protein.
<b>Storage and Stability:</b>	-20°C/1 year