

Antibodies Product Datasheet

Product Name:	CSK mouse mAb
Catalog Number:	N0873
Reactivity:	Human,Rat
Applications:	WB,ICC
Dilutions:	wb 1:1000 icc 1:50
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CSK protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	csk
Protein Name:	
Other Name:	C SRC;C SRC kinase;C src Tyrosine Kinase;C-SRC kinase;c-src tyrosine kinase;Csk A;CSK; CSK_HUMAN;CYTOPLASMIC TYROSINE KINASE;EC 2.7.10.2;MGC112926;MGC117393;MGC154049;P60 Src; Protein tyrosine kinase CYL;Protein-tyrosine kinase CYL;Proto oncogene tyrosine
Human Swiss Prot No.:	P41240
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CSK and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year