

## Antibodies Product Datasheet

<b>Product Name:</b>	Rabbit Anti Human IgG F(ab') <sub>2</sub> - Alexa Fluor 568
<b>Catalog Number:</b>	SA0826
<b>Reactivity:</b>	Human
<b>Conjugate:</b>	Alexa Fluor 568
<b>Applications:</b>	ELISA,IF/ICC,FCM
<b>Dilutions:</b>	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Concentration:</b>	1mg/mL
<b>Other Name:</b>	Rabbit Anti-Human IgG F(ab') <sub>2</sub> - Alexa Fluor 568
<b>Formulation:</b>	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
<b>Storage and Stability:</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing