

Antibodies Product Datasheet

Product Name:	GluR-2 (Phospho-Tyr876) Antibody
Catalog Number:	A9088
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GluR-2 (Phospho-Tyr876)
Concentration:	1 mg/ml
Gene Name:	GRIA2 GLUR2
Protein Name:	GluR-2 (Phospho-Tyr876)
Other Name:	Glutamate receptor 2 (GluR-2) (AMPA-selective glutamate receptor 2) (GluR-B) (GluR-K2) (Glutamate receptor ionotropic, AMPA 2) (GluA2)
Human Swiss Prot No.:	P42262
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of GluR-2 (Phospho-Tyr876) at Human:Y876, Mouse:Y876, Rat:Y876
Storage and Stability:	-20°C/1 year