

Antibodies Product Datasheet

Product Name:	ADRA1A Antibody
Catalog Number:	R12-2456
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	ADRA1A Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADRA1A.
Concentration:	1 mg/ml
Gene Name:	ADRA1A
Other Name:	Alpha-1A adrenergic receptor; Alpha 1A-adrenoceptor; Alpha 1A-adrenoreceptor; Alpha-1C adrenergic receptor; Alpha-adrenergic receptor 1c; ADRA1C; ADA1A
Human Swiss Prot No.:	P35348
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ADRA1A Antibodyantibody detects endogenous levels of ADRA1A.
Storage and Stability:	Stable at -20°C for at least 1 year.