

Antibodies Product Datasheet

Product Name:	ALDH1B1 Antibody
Catalog Number:	C14388
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ALDH1B1. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	ALDH1B1
Protein Name:	Aldehyde dehydrogenase X mitochondrial
Other Name:	ALDH1B1; ALDH5; ALDHX; Aldehyde dehydrogenase X; mitochondrial; Aldehyde dehydrogenase 5; Aldehyde dehydrogenase family 1 member B1
Human Swiss Prot No.:	P30837
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ALDH1B1 Polyclonal Antibody detects endogenous levels of ALDH1B1 protein.
Storage and Stability:	-20°C/1 year