

Antibodies Product Datasheet

Product Name:	p21 (Acetyl-Lys161/K163) Antibody
Catalog Number:	D12216
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human p21 (Acetyl Lys161/K163)
Concentration:	1 mg/ml
Gene Name:	CDKN1A CAP20 CDKN1 CIP1 MDA6 PIC1 SDI1 WAF1
Protein Name:	p21 (Acetyl Lys161/K163)
Other Name:	Cyclin-dependent kinase inhibitor 1 (CDK-interacting protein 1) (Melanoma differentiation-associated protein 6) (MDA-6) (p21)
Human Swiss Prot No.:	P38936
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat p21 (Acetyl Lys161/K163)
Storage and Stability:	-20°C/1 year