

Antibodies Product Datasheet

Product Name:	Caspase-3 Antibody
Catalog Number:	C30154
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CASP3. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	CASP3
Protein Name:	Caspase3
Other Name:	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1
Human Swiss Prot No.:	P42574
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-3 Polyclonal Antibody detects endogenous levels of Caspase-3 protein.
Storage and Stability:	-20°C/1 year