

Antibodies Product Datasheet

Product Name:	SAPAP1 Antibody
Catalog Number:	C10034
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DLGP1. AA range:791-840
Concentration:	1 mg/ml
Gene Name:	DLGAP1
Protein Name:	Disks large-associated protein 1
Other Name:	DLGAP1; DAP1; GKAP; Disks large-associated protein 1; DAP-1; Guanylate kinase-associated protein; hGKAP; PSD-95/SAP90-binding protein 1; SAP90/PSD-95-associated protein 1; SAPAP1
Human Swiss Prot No.:	O14490
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SAPAP1 Polyclonal Antibody detects endogenous levels of SAPAP1 protein.
Storage and Stability:	-20°C/1 year