

Antibodies Product Datasheet

Product Name:	SLC7A6 Antibody
Catalog Number:	C50512
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC7A6
Concentration:	1 mg/ml
Gene Name:	SLC7A6 KIAA0245
Protein Name:	Y+L amino acid transporter 2 (Cationic amino acid transporter, y+ system) (Solute carrier family 7 member 6) (y(+)-type amino acid transporter 2) (Y+LAT2) (y+LAT-2)
Other Name:	
Human Swiss Prot No.:	Q92536
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC7A6 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year