

Antibodies Product Datasheet

Product Name:	MRP-S18C Antibody
Catalog Number:	C14039
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPS18C. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	MRPS18C
Protein Name:	28S ribosomal protein S18c mitochondrial
Other Name:	MRPS18C; CGI-134; 28S ribosomal protein S18c; mitochondrial; MRP-S18-c; Mrps18-c; S18mt-c; 28S ribosomal protein S18-1, mitochondrial; MRP-S18-1
Human Swiss Prot No.:	Q9Y3D5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-S18C Polyclonal Antibody detects endogenous levels of MRP-S18C protein.
Storage and Stability:	-20°C/1 year