

Antibodies Product Datasheet

Product Name:	XRCC5 Antibody
Catalog Number:	R12-3683
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	XRCC5 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human XRCC5.
Concentration:	1 mg/ml
Gene Name:	XRCC5
Other Name:	86 kDa subunit of Ku antigen; ATP-dependent DNA helicase II, 80 kDa subunit; CTC box binding factor 85 kDa subunit; CTC85; CTCBF; DNA-repair protein XRCC5; G22P2; KU86; Lupus Ku autoantigen protein p86; Nuclear factor IV; Thyroid- lupus autoantigen; TLAA;
Human Swiss Prot No.:	P13010
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	XRCC5 Antibody antibody detects endogenous levels of XRCC5.
Storage and Stability:	Stable at -20°C for at least 1 year.