
Antibodies Product Datasheet

Product Name:	eEF2K Antibody
Catalog Number:	R12-2732
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	eEF2K Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human eEF2K.
Concentration:	1 mg/ml
Gene Name:	EEF2K
Other Name:	Calcium/calmodulin-dependent eukaryotic elongation factor-2 kinase; EF2K; Elongation factor 2 kinase; eEF-2 kinase; eEF-2K; kinase eEF2K
Human Swiss Prot No.:	O00418
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	eEF2K Antibody antibody detects endogenous levels of eEF2K.
Storage and Stability:	Stable at -20°C for at least 1 year.