

Antibodies Product Datasheet

Product Name:	CYB5R1 Antibody
Catalog Number:	C16851
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CYB5R1. AA range:9-58
Concentration:	1 mg/ml
Gene Name:	CYB5R1
Protein Name:	NADH-cytochrome b5 reductase 1
Other Name:	CYB5R1; NQO3A2; NADH-cytochrome b5 reductase 1; b5R.1; Humb5R2; NAD(P)H:quinone oxidoreductase type 3 polypeptide A2
Human Swiss Prot No.:	Q9UHQ9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYB5R1 Polyclonal Antibody detects endogenous levels of CYB5R1 protein.
Storage and Stability:	-20°C/1 year