

Antibodies Product Datasheet

Product Name:	Caspase-8 Monoclonal Antibody(2G12)
Catalog Number:	N1510
Reactivity:	Human,Mouse,Rat
Applications:	WB,IF/ICC,IHC-p
Dilutions:	WB: 1:1000-2000 IHC:1:200-500 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of Caspase-8
Concentration:	
Gene Name:	CASP8
Protein Name:	Caspase8
Other Name:	CASP8; MCH5; Caspase-8; CASP-8; Apoptotic cysteine protease; Apoptotic protease Mch-5; CAP4; FADD-homologous ICE/ced-3-like protease; FADD-like ICE; FLICE; ICE-like apoptotic protease 5; MORT1-associated ced-3 homolog; MACH
Human Swiss Prot No.:	Q14790
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Caspase-8 protein.
Storage and Stability:	-20°C/1 year