

Antibodies Product Datasheet

Product Name:	MMP13 (Cleaved-Tyr104) Antibody
Catalog Number:	L0465
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MMP13 (Cleaved-Tyr104)
Concentration:	1 mg/ml
Gene Name:	MMP13
Protein Name:	MMP13 (Cleaved-Tyr104)
Other Name:	Collagenase 3 (EC 3.4.24.-) (Matrix metalloproteinase-13) (MMP-13)
Human Swiss Prot No.:	P45452
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human MMP13 (Cleaved-Tyr104, protein was cleaved amino acid sequence between 103-104)
Storage and Stability:	-20°C/1 year