

Antibodies Product Datasheet

Product Name:	HIF-3 α Antibody
Catalog Number:	R12-2181
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HIF-3 α . AA range:305-354
Concentration:	1 mg/ml
Gene Name:	HIF3A
Protein Name:	Hypoxia-inducible factor 3- α
Other Name:	Hypoxia-inducible factor 3- α ; HIF-3- α ; HIF3- α ; Basic-helix-loop-helix-PAS protein MOP7; Class E basic helix-loop-helix protein 17; bHLHe17; HIF3- α -1; Inhibitory PAS domain protein; IPAS; Member of PAS protein 7; PAS domain-c
Human Swiss Prot No.:	Q66K72
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HIF-3 α Polyclonal Antibody detects endogenous levels of HIF-3 α protein.
Storage and Stability:	-20°C/1 year