

Antibodies Product Datasheet

Product Name:	Cystatin SA Antibody
Catalog Number:	C15282
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CST2. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	CST2
Protein Name:	Cystatin-SA
Other Name:	CST2; Cystatin-SA; Cystatin-2; Cystatin-S5
Human Swiss Prot No.:	P09228
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cystatin SA Polyclonal Antibody detects endogenous levels of Cystatin SA protein.
Storage and Stability:	-20°C/1 year