

Antibodies Product Datasheet

Product Name:	Caspase-7 Antibody
Catalog Number:	L0164
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	WB 1:500-2000, IF 1:50-300, IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Caspase-7 . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	CASP7
Protein Name:	Caspase7
Other Name:	CASP7; MCH3; Caspase-7; CASP-7; Apoptotic protease Mch-3; CMH-1; ICE-like apoptotic protease 3; ICE-LAP3
Human Swiss Prot No.:	P55210
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-7 Polyclonal Antibody detects endogenous levels of Caspase-7 protein.
Storage and Stability:	-20°C/1 year