

## Antibodies Product Datasheet

<b>Product Name:</b>	SPANXB1 Antibody
<b>Catalog Number:</b>	C50055
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit, IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human SPANXB1
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SPANXB1 SPANXB; SPANXB2; SPANXF1
<b>Protein Name:</b>	Sperm protein associated with the nucleus on the X chromosome B/F (Cancer/testis antigen 11.2) (CT11.2) (Nuclear-associated protein SPAN-Xb) (SPANX-B) (Nuclear-associated protein SPAN-Xf) (SPANX-F) (SPANX family member B/F)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q9NS25
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of SPANXB1 at Human
<b>Storage and Stability:</b>	Store at -20° C/1 year