

Antibodies Product Datasheet

Product Name:	14-3-3 ζ Antibody
Catalog Number:	B0001
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IP, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human 14-3-3 zeta. AA range:24-73
Concentration:	1 mg/ml
Gene Name:	YWHAZ
Protein Name:	14-3-3 protein zeta/delta
Other Name:	YWHAZ; 14-3-3 protein zeta/delta; Protein kinase C inhibitor protein 1; KCIP-1
Human Swiss Prot No.:	P63104
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	14-3-3 ζ Polyclonal Antibody detects endogenous levels of 14-3-3 ζ protein.
Storage and Stability:	-20°C/1 year