
Antibodies Product Datasheet

Product Name:	Bcl-2 (Phospho-Thr69) Antibody
Catalog Number:	A0774
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BCL-2 around the phosphorylation site of Thr69. AA range:46-95
Concentration:	1 mg/ml
Gene Name:	BCL2
Protein Name:	Apoptosis regulator Bcl-2
Other Name:	BCL2; Apoptosis regulator Bcl-2
Human Swiss Prot No.:	P10415
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Bcl-2 (T69) Polyclonal Antibody detects endogenous levels of Bcl-2 protein only when phosphorylated at T69.
Storage and Stability:	-20°C/1 year