
Antibodies Product Datasheet

Product Name:	eEF2K Antibody
Catalog Number:	B8340
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human eEF2K. AA range:336-385
Concentration:	1 mg/ml
Gene Name:	EEF2K
Protein Name:	Eukaryotic elongation factor 2 kinase
Other Name:	EEF2K; Eukaryotic elongation factor 2 kinase; eEF-2 kinase; eEF-2K; Calcium/calmodulin-dependent eukaryotic elongation factor 2 kinase
Human Swiss Prot No.:	O00418
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eEF2K Polyclonal Antibody detects endogenous levels of eEF2K protein.
Storage and Stability:	-20°C/1 year