

Antibodies Product Datasheet

Product Name:	C-RAF Antibody
Catalog Number:	R12-2652
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	C-RAF Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C-RAF.
Concentration:	1 mg/ml
Gene Name:	RAF1
Other Name:	C-RAF; C-Raf; CRAF; kinase Raf1; Protein kinase raf 1; RAF proto-oncogene serine/threonine-protein kinase; Raf proto-oncogene serine/threonine-protein kinase; RAF-1; Raf-1; RAF1; Similar to murine leukemia viral (V-raf-1) oncogene homolog 1
Human Swiss Prot No.:	P04049
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	C-RAF Antibody antibody detects endogenous levels of C-RAF.
Storage and Stability:	Stable at -20°C for at least 1 year.