

Antibodies Product Datasheet

Product Name:	TAF5L Antibody
Catalog Number:	C10185
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAF5L. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	TAF5L
Protein Name:	TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L
Other Name:	TAF5L; PAF65B; TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L; PCAF-associated factor 65 beta; PAF65-beta
Human Swiss Prot No.:	O75529
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TAF5L Polyclonal Antibody detects endogenous levels of TAF5L protein.
Storage and Stability:	-20°C/1 year