

Antibodies Product Datasheet

Product Name:	UCH-L1 Antibody
Catalog Number:	C0302
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PGP9.5. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	UCHL1
Protein Name:	Ubiquitin carboxyl-terminal hydrolase isozyme L1
Other Name:	UCHL1; Ubiquitin carboxyl-terminal hydrolase isozyme L1; UCH-L1; Neuron cytoplasmic protein 9.5; PGP 9.5; PGP9.5; Ubiquitin thioesterase L1
Human Swiss Prot No.:	P09936
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	UCH-L1 Polyclonal Antibody detects endogenous levels of UCH-L1 protein.
Storage and Stability:	-20°C/1 year