

Antibodies Product Datasheet

Product Name:	IP3KA Antibody
Catalog Number:	C10462
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IP3KA. AA range:396-445
Concentration:	1 mg/ml
Gene Name:	ІТРКА
Protein Name:	Inositol-trisphosphate 3-kinase A
Other Name:	ITPKA; Inositol-trisphosphate 3-kinase A; Inositol 1; 4,5-trisphosphate 3- kinase A; IP3 3-kinase A; IP3K A; InsP 3-kinase A
Human Swiss Prot No.:	P23677
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IP3KA Polyclonal Antibody detects endogenous levels of IP3KA protein.
Storage and Stability:	-20°C/1 year