

Antibodies Product Datasheet

Product Name:	VEGF-B Antibody
Catalog Number:	C30700
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human VEGFB. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	VEGFB
Protein Name:	Vascular endothelial growth factor B
Other Name:	VEGFB; VRF; Vascular endothelial growth factor B; VEGF-B; VEGF-related factor; VRF
Human Swiss Prot No.:	P49765
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VEGF-B Polyclonal Antibody detects endogenous levels of VEGF-B protein.
Storage and Stability:	-20°C/1 year