

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | AKAP 10 Antibody |
| Catalog Number: | C10110 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human AKAP10. AA range:10-59 |
| Concentration: | 1 mg/ml |
| Gene Name: | AKAP10 |
| Protein Name: | A-kinase anchor protein 10 mitochondrial |
| Other Name: | AKAP10; A-kinase anchor protein 10; mitochondrial; AKAP-10; Dual specificity A kinase-anchoring protein 2; D-AKAP-2; Protein kinase A-anchoring protein 10; PRKA10 |
| Human Swiss Prot No.: | O43572 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | AKAP 10 Polyclonal Antibody detects endogenous levels of AKAP 10 protein. |
| Storage and Stability: | -20°C/1 year |