

## Antibodies Product Datasheet

<b>Product Name:</b>	β III tubulin Monoclonal Antibody(Mix)
<b>Catalog Number:</b>	N1430
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB: 1:5000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of β III tubulin
<b>Concentration:</b>	
<b>Gene Name:</b>	TUBB3
<b>Protein Name:</b>	Tubulin beta-3 chain
<b>Other Name:</b>	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
<b>Human Swiss Prot No.:</b>	Q13509
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous β III tubulin protein.
<b>Storage and Stability:</b>	-20°C/1 year