

Antibodies Product Datasheet

Product Name:	SIPA1 Antibody
Catalog Number:	C51192
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal SIPA1
Concentration:	1 mg/ml
Gene Name:	SIPA1 SPA1
Protein Name:	Signal-induced proliferation-associated protein 1 (Sipa-1) (GTPase-activating protein Spa-1) (p130 SPA-1)
Other Name:	Signal-induced proliferation-associated protein 1 (Sipa-1) (GTPase-activating protein Spa-1) (p130 SPA-1)
Human Swiss Prot No.:	Q96FS4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SIPA1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year