

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | SMC2 Antibody |
| Catalog Number: | C51061 |
| Reactivity: | Human, Mouse |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit, IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human SMC2 |
| Concentration: | 1 mg/ml |
| Gene Name: | SMC2 CAPE SMC2L1 PRO0324 |
| Protein Name: | Structural maintenance of chromosomes protein 2 (SMC protein 2) (SMC-2) (Chromosome-associated protein E) (hCAP-E) (XCAP-E homolog) |
| Other Name: | |
| Human Swiss Prot No.: | O95347 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of SMC2 at Human, Mouse |
| Storage and Stability: | Store at -20° C/1 year |