

Antibodies Product Datasheet

Product Name:	Anti-Human p63 (ABT-P63) Mouse Monoclonal Antibody
Catalog Number:	V0107
Clone ID:	ABT-P63
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	TP63
Gene Name:	Q9H3D4
Other Name:	AIS, Amplified in squamous cell carcinoma, B(p51A), B(p51B), Chronic ulcerative stomatitis protein, CUSP, DN p63 alpha 1, DNp63, EEC3, id:ibd3516, Keratinocyte transcription factor, Keratinocyte transcription factor KET, KET, LMS, MGC115972, MGC192897, NBP, OFC8, OTTHUMP00000209732, OTTHUMP00000209733, OTTHUMP00000209734, OTTHUMP00000209735, OTTHUMP00000209737, OTTHUMP00000209738, OTTHUMP00000209739, OTTHUMP00000209740, OTTHUMP00000209741, OTTHUMP00000209742, OTTHUMP00000209743, OTTHUMP00000209744, p40, p51, P51/P63, p53-related protein p63, p53CP, p63, P63_HUMAN, p73H, p73L, RHS, SHFM4, TAp63alpha, TP53CP, TP53L, TP63, TP73L, Transformation related protein 63, Transformation-related protein 63, Trp53rp1, Trp63, Tumor protein 63, Tumor protein p53-competing protein, Tumor protein p53-like, Tumor protein p63, Tumor protein p63 deltaN isoform delta, Tumor protein p73, Tumor protein p73-like
Isotype	IgG1, Kappa
Localization:	Tonsil
Positive Tissue Control:	Nuclear
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human p63 protein, including ΔNp63 and TAp63. In western blotting of A431, Hela cell lysates, the antibody can label a 77 kDa band corresponding to p63, while there is no band detected in-

V0107 | Anti-Human p63 (ABT-P63) Mouse Monoclonal Antibody

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022

p63 knock out Hela cell lys

Storage and Stability: -20°C/1 year