

Antibodies Product Datasheet

Product Name:	Stat1 (Acetyl-Lys410/K413) Antibody
Catalog Number:	D12320
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Stat1 (Acetyl Lys410/K413)
Concentration:	1 mg/ml
Gene Name:	STAT1
Protein Name:	Stat1 (Acetyl Lys410/K413)
Other Name:	Signal transducer and activator of transcription 1-alpha/beta (Transcription factor ISGF-3 components p91/p84)
Human Swiss Prot No.:	P42224
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Stat1 (Acetyl Lys410/K413)
Storage and Stability:	-20°C/1 year