
Antibodies Product Datasheet

Product Name:	ApoO Monoclonal Antibody
Catalog Number:	N0040
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF (paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of ApoO expressed in E. Coli.
Concentration:	
Gene Name:	APOO
Protein Name:	Apolipoprotein O
Other Name:	APOO; FAM121B; My025; Apolipoprotein O; Protein FAM121B
Human Swiss Prot No.:	Q9BUR5
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	ApoO Monoclonal Antibody detects endogenous levels of ApoO protein.
Storage and Stability:	-20°C/1 year