

Antibodies Product Datasheet

Product Name:	ADAM32 Antibody
Catalog Number:	C14296
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAM32 . at AA range: 460-540
Concentration:	1 mg/ml
Gene Name:	ADAM32
Protein Name:	Disintegrin and metalloproteinase domain-containing protein 32
Other Name:	ADAM32; Disintegrin and metalloproteinase domain-containing protein 32; ADAM 32
Human Swiss Prot No.:	Q8TC27
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAM32 Polyclonal Antibody detects endogenous levels of ADAM32 protein.
Storage and Stability:	-20°C/1 year