

Antibodies Product Datasheet

Product Name:	YIPF5 Antibody
Catalog Number:	C50326
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human YIPF5
Concentration:	1 mg/ml
Gene Name:	YIPF5 FINGER5 YIP1A PP12723 SB140 UNQ3123/PRO10275
Protein Name:	Protein YIPF5 (Five-pass transmembrane protein localizing in the Golgi apparatus and the endoplasmic reticulum 5) (Smooth muscle cell-associated protein 5) (SMAP-5) (YIP1 family member 5) (YPT-interacting protein 1 A)
Other Name:	
Human Swiss Prot No.:	Q969M3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of YIPF5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year