

Assay Biotechnology Company

47787 Fremont Blvd.
Fremont, CA 94538
Phone: (408) 747 0185
Fax: (408 747 0145
www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: ARK-1 Antibody

Catalog Number: B1132

Reactivity: Human, Mouse, Rat, Monkey

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in

other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human AurA. AA range:311-360

Concentration: 1 mg/ml

Gene Name: AURKA

Protein Name: Aurora kinase A

Other Name: AURKA; AIK; AIRK1; ARK1; AURA; AYK1; BTAK; IAK1; STK15; STK6;

Aurora kinase A; Aurora 2; Aurora/IPL1-related kinase 1; ARK-1; Aurora-

related kinase 1; hARK1; Breast tumor-amplified kinase;

Serine/threonine-protein kinase 15; Serine/threonin

Human Swiss Prot No.: 014965

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: ARK-1 Polyclonal Antibody detects endogenous levels of ARK-1 protein.

Storage and Stability: -20°C/1 year