

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Nup93 Antibody |
| Catalog Number: | C17097 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human NUP93. AA range:221-270 |
| Concentration: | 1 mg/ml |
| Gene Name: | NUP93 |
| Protein Name: | Nuclear pore complex protein Nup93 |
| Other Name: | NUP93; KIAA0095; Nuclear pore complex protein Nup93; 93 kDa nucleoporin; Nucleoporin Nup93 |
| Human Swiss Prot No.: | Q8N1F7 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Nup93 Polyclonal Antibody detects endogenous levels of Nup93 protein. |
| Storage and Stability: | -20°C/1 year |