

Antibodies Product Datasheet

Product Name:	PDGFRb (Phospho-Tyr857) Antibody
Catalog Number:	A8167
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PDGFRb (Phospho-Tyr857)
Concentration:	1 mg/ml
Gene Name:	PDGFRB PDGFR PDGFR1
Protein Name:	PDGFRb (Phospho-Tyr857)
Other Name:	Platelet-derived growth factor receptor beta (PDGF-R-beta) (PDGFR-beta) (EC 2.7.10.1) (Beta platelet-derived growth factor receptor) (Beta-type platelet-derived growth factor receptor) (CD140 antigen-like family member B) (Platelet-derived growth factor r
Human Swiss Prot No.:	P09619
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of PDGFRb (Phospho-Tyr857) at Human:Y857, Mouse:Y856, Rat:Y856
Storage and Stability:	-20°C/1 year