

Antibodies Product Datasheet

Product Name:	ERF Antibody
Catalog Number:	C10666
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ERF. AA range:471-520
Concentration:	1 mg/ml
Gene Name:	ERF
Protein Name:	ETS domain-containing transcription factor ERF
Other Name:	ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2
Human Swiss Prot No.:	P50548
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERF Polyclonal Antibody detects endogenous levels of ERF protein.
Storage and Stability:	-20°C/1 year