
Antibodies Product Datasheet

Product Name:	GOLT1A Antibody
Catalog Number:	C19496
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GOLT1A. AA range:40-89
Concentration:	1 mg/ml
Gene Name:	GOLT1A
Protein Name:	Vesicle transport protein GOT1A
Other Name:	GOLT1A; GOT1B; Vesicle transport protein GOT1A; Golgi transport 1 homolog A; hGOT1b
Human Swiss Prot No.:	Q6ZVE7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GOLT1A Polyclonal Antibody detects endogenous levels of GOLT1A protein.
Storage and Stability:	-20°C/1 year